



# DR400 patch style recorder

## Revision Chart

	Rev A	Rev B
		
Started shipping	July 2019	December 2020
Charging and uploading	USB into back of recorder	Docking station
Firmware	5.03 – 5.07	5.21
PC Patch version	Can use any version	Must use version 1.08
Button on front says	EVENT	PRESS
When plugged into PC, LED light ...	Indicates PCPatch responses	Indicates charging status (below)
When plugged into power, LED light...	Flashes orange when charging. Flashes green when fully charged.	



### NorthEast Monitoring, Inc.

141 Parker Street, Suite 200  
 Maynard, MA 01754 USA  
[www.nemon.com](http://www.nemon.com)

phone: [+1] 978-461-3992  
 fax: [+1] 978-461-5991  
 toll free USA: 866-346-5837  
 email: [info@nemon.com](mailto:info@nemon.com)

*advancing Holter technology*